第6版
Sixth Edition

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Laboratory medicine

Andrology Laboratory Director

ANOVA, Clinic for Endocrinology, Karolinska University Hospital and Department of Medicine, Huddinge, Karolinska Institutet, Stockholm, Sweden

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Professor Oleg Apolikhin, MD, PhD

Andrology

Professor and Director, Research Institute of Urology and Interventional Radiology named after N.A. Lopatkin-branch of the National Medical Research Center of Radiology of the Ministry of Health of the Russian Federation
Moscow, Russian Federation

Professor Elisabetta Baldi, Doctor in Biology

Laboratory medicine

Associate Professor of Clinical Pathology, Department of Experimental and Clinical Medicine, University of Florence

Responsible of the semen cryopreservation centre, Azienda Ospedaliero-Universitaria Careggi (AOUC)

Florence, Italy

Professor Christopher L.R. Barratt PhD, DSc, FRSE

Head of Division Systems Medicine, Professor of Reproductive Medicine

School of Medicine, University of Dundee

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Professor Mario Philip R. Festin, MD, MSc, MHPed

Obstetrics-gynaecology, clinical epidemiology, health professions education, family planning
Professor, Department of Obstetrics and Gynecology and Department of Clinical Epidemiology,
College of Medicine, University of the Philippines Manila
Convenor, Philippine Study Group on Reproductive Health, National Institutes of Health,
University of the Philippines
Manila, Philippines
WHO Secretariat Consultant

Dr. Jackson C. Kirkman-Brown, MBE, PhD

Reproductive science, andrology, health care science, fertility
University Reader and Director, Centre for Human Reproductive Science (ChRS), College of
Medical and Dental Sciences, University of Birmingham
HFEA Person Responsible and Science Lead, Birmingham Women's Fertility Centre, Birmingham
Women's and Children's NHS Foundation Trust
Birmingham, United Kingdom

Professor Dolores J. Lamb, MS, PhD, HCLD (ABB)

Andrology/urology, clinical chemistry (National Registry of Clinical Chemistry)
High Complexity Laboratory Director (HCLD, American Board of Bioanalysts)
Dow Professor of Urology
Vice-Chairman for Research (Urology)
James Buchanan Brady Foundation Department of Urology, Director, Center for Reproductive
Genomics, Englander Institute for Precision Medicine
Weill Cornell Medical College
New York, United States of America
American Society for Reproductive Medicine

Dr. Michael Mbizvo, PhD

Andrology and reproductive epidemiology
Professor, Reproductive Health Sciences, University of Zimbabwe
Country Director/Senior Associate, Population Council
Lusaka, Zambia

Prof. Dr. Stefan Schlatt, PhD

Reproductive biology, andrology
University Professor and Director, Centre for Reproductive Medicine and Andrology, University
of Münster, Germany
WHO Collaborating Centre for Research in Male Reproduction

Dr. Igor Toskin, MD, PhD, DSc

Public health, STI, HIV/AIDS, urology, dermatology, epidemiology, sexual and reproductive health
Scientist, Department of Sexual and Reproductive Health and Research
World Health Organization
Geneva, Switzerland
WHO Secretariat

Dr. Christina Wang, MD

Endocrinology, andrology (male reproduction and contraception)
Professor of Medicine, David Geffen School of Medicine, UCLA
Division of Endocrinology, Department of Medicine, Harbor-UCLA Medical Center and the
Lundquist Institute
Torrance California, United States of America

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Karel Blondeel, MSc, MA

HIV/STI prevention and sexual health
Doctoral student, Department of Medicine and Health Sciences, Ghent University
Ghent, Belgium
WHO Secretariat Consultant

MUDr. Petr Houska

Laboratory medicine
Doctoral student, ANOVA, Clinic for Endocrinology, Karolinska University Hospital and
Department of Medicine, Huddinge, Karolinska Institutet
Stockholm, Sweden

Dr. Stepan Krasnyak, MD

Andrology, urology
Associate Scientist, Department of Andrology and Human Reproduction, Research Institute of
Urology and Interventional Radiology named after N.A. Lopatkin-branch of the National Medical
Research Center of Radiology of the Ministry of Health of the Russian Federation Moscow,
Russian Federation

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1. Belsey M, Eliasson R, Gallegos AJ, Moghissi KS, Paulsen CA, Prasad AMN. Laboratory Manual for the Examination of Human Semen and Semen-Cervical Mucus Interaction. Singapore: Press Concern; 1980.
2. World Health Organization. WHO laboratory manual for the examination of human semen and semen-cervical mucus interactions. 2nd ed. Cambridge, UK: Cambridge University Press; 1987. 67 p.
3. World Health Organization. WHO laboratory manual for the examination of human semen and sperm-cervical mucus interactions. 3rd ed. Cambridge, UK: Cambridge University Press; 1992. 107 p.
4. World Health Organization. WHO laboratory manual for the examination of human semen and sperm-cervical mucus interactions. 4th ed. Cambridge, UK: Cambridge University Press; 1999. 128 p.
5. Campbell MJ, Lotti F, Baldi E, Schlatt S, Festin MP, Björndahl L et al. Distribution of Semen Examination Results 2020 - a follow up of data collated for the WHO semen analysis manual 2010. Andrology. 2021.
6. Björndahl L. What is normal semen quality? On the use and abuse of reference limits for the interpretation of semen analysis results. Hum Fertil (Camb). 2011;14(3):179-86.
7. Jedrzejczak P, Taszarek-Hauke G, Hauke J, Pawelczyk L, Duleba AJ. Prediction of spontaneous conception based on semen parameters. Int J Androl. 2008;31(5):499-507.
8. Guzick DS, Overstreet JW, Factor-Litvak P, Brazil CK, Nakajima ST, Coutifaris C et al. Sperm morphology, motility, and concentration in fertile and infertile men. N Engl J Med. 2001;345(19):1388-93.
9. Barratt CLR, Björndahl L, De Jonge CJ, Lamb DJ, Osorio Martini F, McLachlan R et al. The diagnosis of male infertility: an analysis of the evidence to support the development of global WHO guidance-challenges and future research opportunities. Hum Reprod Update. 2017;23(6):660-80.
10. Wang C, Swerdloff RS. Limitations of semen analysis as a test of male fertility and anticipated needs from newer tests. Fertil Steril. 2014;102(6):1502-7.

11. Nieschlag E, Pock T, Hellenkemper B. External quality of semen analysis reveals low compliance with the WHO guidelines. *J Reproduktionsmed Endokrinol_Online*. 2017;14(6):306-10.
12. Barratt CL, Björndahl L, Menkveld R, Mortimer D. ESHRE special interest group for andrology basic semen analysis course: a continued focus on accuracy, quality, efficiency and clinical relevance. *Hum Reprod*. 2011;26(12):3207-12.
13. Björndahl L, Kvist U. Sequence of ejaculation affects the spermatozoon as a carrier and its message. *Reprod Biomed Online*. 2003;7(4):440-8.
14. Weiske WH. Pregnancy caused by sperm from vasa efferentia. *Fertil Steril*. 1994;62(3):642-3.
15. Sobrero AJ, Macleod J. The immediate postcoital test. *Fertil Steril*. 1962;13:184-9.
16. Cooper TG, Barfield JP, Yeung CH. Changes in osmolality during liquefaction of human semen. *Int J Androl*. 2005;28(1):58-60.
17. Holmes E, Björndahl L, Kvist U. Post-ejaculatory increase in human semen osmolality in vitro. *Andrologia*. 2019:e13311.
18. Holmes E, Björndahl L, Kvist U. Possible factors influencing post-ejaculatory changes of the osmolality of human semen in vitro. *Andrologia*. 2019:e13443.
19. Zavos PM, Goodpasture JC. Clinical improvements of specific seminal deficiencies via intercourse with a seminal collection device versus masturbation. *Fertil Steril*. 1989;51(1):190-3.
20. Pound N, Javed MH, Ruberto C, Shaikh MA, Del Valle AP. Duration of sexual arousal predicts semen parameters for masturbatory ejaculates. *Physiology & behavior*. 2002;76(4-5):685-9.
21. Eliasson R. Basic semen analysis. In: Matson P, editor. *Current Topics in Andrology*. Perth: Ladybrook Publishing; 2003. p. 35-89.
22. Amann RP. Evaluating spermatogenesis using semen: the biology of emission tells why reporting total sperm per sample is important, and why reporting only number of sperm per milliliter is irrational. *J Androl*. 2009;30(6):623-5.
23. Ng KK, Donat R, Chan L, Lalak A, Di Pierro I, Handelsman DJ. Sperm output of older men. *Hum Reprod*. 2004;19(8):1811-5.
24. Cooper TG, Keck C, Oberdieck U, Nieschlag E. Effects of multiple ejaculations after extended periods of sexual abstinence on total, motile and normal sperm numbers, as well as accessory gland secretions, from healthy normal and oligozoospermic men. *Hum Reprod*. 1993;8(8):1251-8.

25. Fernandez JL, Muriel L, Goyanes V, Segrelles E, Gosalvez J, Enciso M et al. Simple determination of human sperm DNA fragmentation with an improved sperm chromatin dispersion test. *Fertil Steril*. 2005;84(4):833-42.
26. Tyler JP, Crockett NG, Driscoll GL. Studies of human seminal parameters with frequent ejaculation. I. Clinical characteristics. *Clin Reprod Fertil*. 1982;1(4):273-85.
27. Tyler JP, Crockett NG, Driscoll GL. Studies of human seminal parameters with frequent ejaculation II. Spermatozoal vitality and storage. *Clin Reprod Fertil*. 1982;1(4):287-93.
28. De Jonge C, LaFromboise M, Bosmans E, Ombelet W, Cox A, Nijs M. Influence of the abstinence period on human sperm quality. *Fertil Steril*. 2004;82(1):57-65.
29. Correa-Perez JR. Intrasubject abstinence and quality of sperm. *Fertil Steril*. 2005;83(5):1588-9; author reply 9-90.
30. Amann RP, Howards SS. Daily spermatozoal production and epididymal spermatozoal reserves of the human male. *J Urol*. 1980;124(2):211-5.
31. Amann RP. Considerations in evaluating human spermatogenesis on the basis of total sperm per ejaculate. *J Androl*. 2009;30(6):626-41.
32. Handelsman DJ, Conway AJ, Boylan LM, Turtle JR. Testicular function in potential sperm donors: normal ranges and the effects of smoking and varicocele. *Int J Androl*. 1984;7(5):369-82.
33. Behre HM, Bergmann M, Simoni M, Tuttelmann F. Primary Testicular Failure. In: De Groot LJ, Chrousos G, Dungan K, Feingold KR, Grossman A, Hershman JM, et al., editors. *Endotext*. South Dartmouth (MA)2000.
34. Andersen AG, Jensen TK, Carlsen E, Jorgensen N, Andersson AM, Krarup T et al. High frequency of sub-optimal semen quality in an unselected population of young men. *Hum Reprod*. 2000;15(2):366-72.
35. Holstein AF, Schulze W, Davidoff M. Understanding spermatogenesis is a prerequisite for treatment. *Reprod Biol Endocrinol*. 2003;1:107.
36. Santi D, Corona G. Primary and Secondary Hypogonadism. In: Simoni M, Huhtaniemi I, editors. *Endocrinology of the Testis and Male Reproduction*. 2. London: Springer International Publishing; 2017. p. 687-747.
37. Medina P, Segarra G, Ballester R, Chuan P, Domenech C, Vila JM et al. Effects of antidepressants in adrenergic neurotransmission of human vas deferens. *Urology*. 2000;55(4):592-7.
38. Baker HW, Kovacs GT. Spontaneous improvement in semen quality: regression towards the mean. *Int J Androl*. 1985;8(6):421-6.

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39. Alvarez C, Castilla JA, Martinez L, Ramirez JP, Vergara F, Gaforio JJ. Biological variation of seminal parameters in healthy subjects. *Hum Reprod.* 2003;18(10):2082-8.
40. Poland ML, Moghissi KS, Giblin PT, Ager JW, Olson JM. Variation of semen measures within normal men. *Fertil Steril.* 1985;44(3):396-400.
41. Berman NG, Wang C, Paulsen CA. Methodological issues in the analysis of human sperm concentration data. *J Androl.* 1996;17(1):68-73.
42. Carlsen E, Petersen JH, Andersson AM, Skakkebaek NE. Effects of ejaculatory frequency and season on variations in semen quality. *Fertil Steril.* 2004;82(2):358-66.
43. Castilla JA, Alvarez C, Aguilar J, Gonzalez-Varea C, Gonzalvo MC, Martinez L. Influence of analytical and biological variation on the clinical interpretation of seminal parameters. *Hum Reprod.* 2006;21(4):847-51.
44. Keel BA. Within- and between-subject variation in semen parameters in infertile men and normal semen donors. *Fertil Steril.* 2006;85(1):128-34.
45. Mortimer D. *Practical Laboratory Andrology.* Oxford: Oxford University Press; 1994. 393 p.
46. Mortimer D. Laboratory standards in routine clinical andrology. *Reproductive Medicine Review.* 1994;3:97-111.
47. Jones DM, Kovacs GT, Harrison L, Jennings MG, Baker HW. Immobilization of sperm by condoms and their components. *Clin Reprod Fertil.* 1986;4(6):367-72.
48. World Health Organization. *Laboratory biosafety manual*, 3rd edition. Geneva 2004.
49. Huggins C, Scott WW, Heinen JH. Chemical composition of human semen and of the secretions of the prostate and seminal vesicles. *Am J Physiol.* 1942;136(3):467-73.
50. Brazil C, Swan SH, Drobnis EZ, Liu F, Wang C, Redmon JB et al. Standardized methods for semen evaluation in a multicenter research study. *J Androl.* 2004;25(4):635-44.
51. Cooper TG, Brazil C, Swan SH, Overstreet JW. Ejaculate volume is seriously underestimated when semen is pipetted or decanted into cylinders from the collection vessel. *J Androl.* 2007;28(1):1-4.
52. Rose NR, Hjort T, Rümke P, Harper MJK, Vyazov O. Techniques for detection of iso- and auto-antibodies to human spermatozoa. *Clin Exp Immunol.* 1976;23(2):175-99.

53. Johannisson E, Campana A, Luthi R, de Agostini A. Evaluation of 'round cells' in semen analysis: a comparative study. *Hum Reprod Update*. 2000;6(4):404-12.
54. Freund M. Standards for the rating of human sperm morphology. A cooperative study. *Int J Fertil*. 1966;11(1):97-180.
55. Jue JS, Ramasamy R. Significance of positive semen culture in relation to male infertility and the assisted reproductive technology process. *Transl Androl Urol*. 2017;6(5):916-22.
56. Micillo A, Vassallo MR, Cordeschi G, D'Andrea S, Necozione S, Francavilla F et al. Semen leukocytes and oxidative-dependent DNA damage of spermatozoa in male partners of subfertile couples with no symptoms of genital tract infection. *Andrology*. 2016;4(5):808-15.
57. Palermo GD, Neri QV, Cozzubbo T, Cheung S, Pereira N, Rosenwaks Z. Shedding Light on the Nature of Seminal Round Cells. *PLoS One*. 2016;11(3):e0151640.
58. Jouannet P, Ducot B, Feneux D, Spira A. Male factors and the likelihood of pregnancy in infertile couples. I. Study of sperm characteristics. *Int J Androl*. 1988;11(5):379-94.
59. Zinaman MJ, Brown CC, Selevan SG, Clegg ED. Semen quality and human fertility: a prospective study with healthy couples. *J Androl*. 2000;21(1):145-53.
60. Larsen L, Scheike T, Jensen TK, Bonde JP, Ernst E, Hjollund NH et al. Computer-assisted semen analysis parameters as predictors for fertility of men from the general population. The Danish First Pregnancy Planner Study Team. *Hum Reprod*. 2000;15(7):1562-7.
61. Aitken RJ, Sutton M, Warner P, Richardson DW. Relationship between the movement characteristics of human spermatozoa and their ability to penetrate cervical mucus and zona-free hamster oocytes. *J Reprod Fertil*. 1985;73(2):441-9.
62. Mortimer D, Pandya IJ, Sawers RS. Relationship between human sperm motility characteristics and sperm penetration into human cervical mucus in vitro. *J Reprod Fertil*. 1986;78(1):93-102.
63. Comhaire FH, Vermeulen L, Hinting A, Schoonjans F. Accuracy of sperm characteristics in predicting the in vitro fertilizing capacity of semen. *J In Vitro Fert Embryo Transf*. 1988;5(6):326-31.
64. Barratt CL, McLeod ID, Dunphy BC, Cooke ID. Prognostic value of two putative sperm function tests: hypo-osmotic swelling and bovine sperm mucus penetration test (Penetrak). *Hum Reprod*. 1992;7(9):1240-4.
65. Irvine DS, Aitken RJ. Predictive value of in-vitro sperm function tests in the context of an AID service. *Hum Reprod*. 1986;1(8):539-45.

66. Bollendorf A, Check JH, Lurie D. Evaluation of the effect of the absence of sperm with rapid and linear progressive motility on subsequent pregnancy rates following intrauterine insemination or in vitro fertilization. *J Androl.* 1996;17(5):550-7.
67. Sifer C, Sasportes T, Barraud V, Poncelet C, Rudant J, Porcher R et al. World Health Organization grade 'a' motility and zona-binding test accurately predict IVF outcome for mild male factor and unexplained infertilities. *Hum Reprod.* 2005;20(10):2769-75.
68. Van den Bergh M, Emiliani S, Biramane J, Vannin AS, Englert Y. A first prospective study of the individual straight line velocity of the spermatozoon and its influences on the fertilization rate after intracytoplasmic sperm injection. *Hum Reprod.* 1998;13(11):3103-7.
69. Björndahl L. The usefulness and significance of assessing rapidly progressive spermatozoa. *Asian J Androl.* 2010;12(1):33-5.
70. Eliasson R. Semen analysis with regard to sperm number, sperm morphology and functional aspects. *Asian J Androl.* 2010;12(1):26-32.
71. Afzelius BA, Eliasson R, Johnsen O, Lindholmer C. Lack of dynein arms in immotile human spermatozoa. *J Cell Biol.* 1975;66(2):225-32.
72. Chemes EH, Rawe YV. Sperm pathology: a step beyond descriptive morphology. Origin, characterization and fertility potential of abnormal sperm phenotypes in infertile men. *Hum Reprod Update.* 2003;9(5):405-28.
73. Wilton LJ, Temple-Smith PD, Baker HW, de Kretser DM. Human male infertility caused by degeneration and death of sperm in the epididymis. *Fertil Steril.* 1988;49(6):1052-8.
74. Correa-Perez JR, Fernandez-Pelegrina R, Aslanis P, Zavos PM. Clinical management of men producing ejaculates characterized by high levels of dead sperm and altered seminal plasma factors consistent with epididymal necropermia. *Fertil Steril.* 2004;81(4):1148-50.
75. Björndahl L, Soderlund I, Johansson S, Mohammadi M, Pourian MR, Kvist U. Why the WHO recommendations for eosin-nigrosin staining techniques for human sperm vitality assessment must change. *J Androl.* 2004;25(5):671-8.
76. Björndahl L, Soderlund I, Kvist U. Evaluation of the one-step eosin-nigrosin staining technique for human sperm vitality assessment. *Hum Reprod.* 2003;18(4):813-6.
77. Mortimer D. A technical note on the assessment of human sperm vitality using eosin-nigrosin staining. *Reprod Biomed Online.* 2020;40(6):851-5.

78. Slama R, Eustache F, Ducot B, Jensen TK, Jorgensen N, Horte A et al. Time to pregnancy and semen parameters: a cross-sectional study among fertile couples from four European cities. *Hum Reprod.* 2002;17(2):503-15.
79. World Health Organization. Task Force on Methods for the Regulation of Male Fertility: Contraceptive efficacy of testosterone-induced azoospermia and oligozoospermia in normal men. *Fertil Steril.* 1996;65(4):821-9.
80. Bonde JP, Ernst E, Jensen TK, Hjollund NH, Kolstad H, Henriksen TB et al. Relation between semen quality and fertility: a population-based study of 430 first-pregnancy planners. *Lancet.* 1998;352(9135):1172-7.
81. Behre HM, Yeung CH, Holstein AF, Weinbauer GF, Gassner C, Nieschlag E. Diagnosis of Male Infertility and Hypogonadism. In: Nieschlag E, Behre HM, editors. *Andrology: Male Reproductive Health and Dysfunction.* second ed. Berlin Heidelberg New York: Springer Verlag; 2000. p. 90-121.
82. MacLeod J, Wang Y. Male fertility potential in terms of semen quality: a review of the past, a study of the present. *Fertil Steril.* 1979;31(2):103-16.
83. Eliasson R. Analysis of semen. In: Behrman SJ, Kistner RW, editors. *Progress in Infertility.* Boston: Little, Brown and Co; 1975. p. 691-713.
84. Björndahl L, Mortimer D, Barratt CLR, Castilla JA, Menkveld R, Kvist U et al. *A Practical Guide to Basic Laboratory Andrology.* Cambridge: Cambridge University Press; 2010.
85. Mortimer D, Shu MA, Tan R. Standardization and quality control of sperm concentration and sperm motility counts in semen analysis. *Human reproduction (Oxford, England).* 1986;1(5):299-303.
86. Kirkman-Brown J, Björndahl L. Evaluation of a disposable plastic Neubauer counting chamber for semen analysis. *Fertil Steril.* 2009;91(2):627-31.
87. Cooper TG, Hellenkemper B, Jonckheere J, Callewaert N, Grootenhuis AJ, Kersemaekers WM et al. Azoospermia: virtual reality or possible to quantify? *J Androl.* 2006;27(4):483-90.
88. Ezech UI, Moore HM. Redefining azoospermia and its implications. *Fertil Steril.* 2001;75(1):213-4.
89. Sharif K. Reclassification of azoospermia: the time has come? *Hum Reprod.* 2000;15(2):237-8.
90. Eliasson R. Analysis of semen. In: Burger HG, De Kretser DM, editors. *The Testis.* New York: Raven Press; 1981. p. 381-99.
91. Lindsay KS, Floyd I, Swan R. Classification of azoospermic samples. *Lancet.* 1995;345(8965):1642.

92. Jaffe TM, Kim ED, Hoekstra TH, Lipshultz LI. Sperm pellet analysis: a technique to detect the presence of sperm in men considered to have azoospermia by routine semen analysis. *J Urol.* 1998;159(5):1548-50.
93. Corea M, Campagnone J, Sigman M. The diagnosis of azoospermia depends on the force of centrifugation. *Fertil Steril.* 2005;83(4):920-2.
94. Fredricsson B, Bjork G. Morphology of postcoital spermatozoa in the cervical secretion and its clinical significance. *Fertil Steril.* 1977;28(8):841-5.
95. Menkveld R, Stander FS, Kotze TJ, Kruger TF, van Zyl JA. The evaluation of morphological characteristics of human spermatozoa according to stricter criteria. *Hum Reprod.* 1990;5(5):586-92.
96. Menkveld R, Franken DR, Kruger TF, Oehninger S, Hodgen GD. Sperm selection capacity of the human zona pellucida. *Mol Reprod Dev.* 1991;30(4):346-52.
97. Liu DY, Baker HW. Morphology of spermatozoa bound to the zona pellucida of human oocytes that failed to fertilize in vitro. *J Reprod Fertil.* 1992;94(1):71-84.
98. Eggert-Kruse W, Schwarz H, Rohr G, Demirakca T, Tilgen W, Runnebaum B. Sperm morphology assessment using strict criteria and male fertility under in-vivo conditions of conception. *Hum Reprod.* 1996;11(1):139-46.
99. Coetzee K, Kruger TF, Lombard CJ. Predictive value of normal sperm morphology: a structured literature review. *Hum Reprod Update.* 1998;4(1):73-82.
100. Toner JP, Mossad H, Grow DR, Morshedi M, Swanson RJ, Oehninger S. Value of sperm morphology assessed by strict criteria for prediction of the outcome of artificial (intrauterine) insemination. *Andrologia.* 1995;27(3):143-8.
101. Van Waart J, Kruger TF, Lombard CJ, Ombelet W. Predictive value of normal sperm morphology in intrauterine insemination (IUI): a structured literature review. *Hum Reprod Update.* 2001;7(5):495-500.
102. Garrett C, Liu DY, Clarke GN, Rushford DD, Baker HW. Automated semen analysis: 'zona pellucida preferred' sperm morphometry and straight-line velocity are related to pregnancy rate in subfertile couples. *Hum Reprod.* 2003;18(8):1643-9.
103. Liu DY, Garrett C, Baker HW. Low proportions of sperm can bind to the zona pellucida of human oocytes. *Hum Reprod.* 2003;18(11):2382-9.
104. Menkveld R, Wong WY, Lombard CJ, Wetzels AM, Thomas CM, Merkus HM et al. Semen parameters, including WHO and strict criteria morphology, in a fertile and subfertile population: an effort towards standardization of in-vivo thresholds. *Hum Reprod.* 2001;16(6):1165-71.

105. van der Merwe FH, Kruger TF, Oehninger SC, Lombard CJ. The use of semen parameters to identify the subfertile male in the general population. *Gynecol Obstet Invest.* 2005;59(2):86-91.
106. Mortimer D. The functional anatomy of the human spermatozoon: relating ultrastructure and function. *Mol Hum Reprod.* 2018;24(12):567-92.
107. Katz DF, Overstreet JW, Samuels SJ, Niswander PW, Bloom TD, Lewis EL. Morphometric analysis of spermatozoa in the assessment of human male fertility. *J Androl.* 1986;7(4):203-10.
108. Abraham-Peskir JV, Chantler E, Uggerhoj E, Fedder J. Response of midpiece vesicles on human sperm to osmotic stress. *Hum Reprod.* 2002;17(2):375-82.
109. Cooper TG, Yeung CH, Fetic S, Sobhani A, Nieschlag E. Cytoplasmic droplets are normal structures of human sperm but are not well preserved by routine procedures for assessing sperm morphology. *Hum Reprod.* 2004;19(10):2283-8.
110. Kruger TF, Acosta AA, Simmons KF, Swanson RJ, Matta JF, Oehninger S. Predictive value of abnormal sperm morphology in in vitro fertilization. *Fertil Steril.* 1988;49(1):112-7.
111. Zollner U, Schleyer M, Steck T. Evaluation of a cut-off value for normal sperm morphology using strict criteria to predict fertilization after conventional in-vitro fertilization and embryo transfer in asthenozoospermia. *Hum Reprod.* 1996;11(10):2155-61.
112. Franken DR, Barendsen R, Kruger TF. A continuous quality control program for strict sperm morphology. *Fertil Steril.* 2000;74(4):721-4.
113. Menkveld R, El-Garem Y, Schill WB, Henkel R. Relationship between human sperm morphology and acrosomal function. *J Assist Reprod Genet.* 2003;20(10):432-8.
114. Kvist U, Björndahl L. Manual on Basic Semen Analysis 2002. Heineman MJ, editor. Oxford: Oxford Academic; 2002 2002-06-15. 24 p.
115. Kruger TF, Menkveld R, Stander FS, Lombard CJ, Van der Merwe JP, van Zyl JA et al. Sperm morphologic features as a prognostic factor in in vitro fertilization. *Fertil Steril.* 1986;46(6):1118-23.
116. Inaba K, Mizuno K. Sperm dysfunction and ciliopathy. *Reprod Med Biol.* 2016;15(2):77-94.
117. Gandini L, Lombardo F, Paoli D, Caponecchia L, Familiari G, Verlengia C et al. Study of apoptotic DNA fragmentation in human spermatozoa. *Hum Reprod.* 2000;15(4):830-9.

118. Lee JD, Kamiguchi Y, Yanagimachi R. Analysis of chromosome constitution of human spermatozoa with normal and aberrant head morphologies after injection into mouse oocytes. *Hum Reprod.* 1996;11(9):1942-6.
119. Dadoune JP, Mayaux MJ, Guihard-Moscato ML. Correlation between defects in chromatin condensation of human spermatozoa stained by aniline blue and semen characteristics. *Andrologia.* 1988;20(3):211-7.
120. Devillard F, Metzler-Guillemain C, Pelletier R, DeRobertis C, Bergues U, Hennebicq S et al. Polyploidy in large-headed sperm: FISH study of three cases. *Hum Reprod.* 2002;17(5):1292-8.
121. Martin RH, Rademaker AW, Greene C, Ko E, Hoang T, Barclay L et al. A comparison of the frequency of sperm chromosome abnormalities in men with mild, moderate, and severe oligozoospermia. *Biol Reprod.* 2003;69(2):535-9.
122. Mortimer D, Menkveld R. Sperm morphology assessment--historical perspectives and current opinions. *J Androl.* 2001;22(2):192-205.
123. Chantler E, Abraham-Peskir JV. Significance of midpiece vesicles and functional integrity of the membranes of human spermatozoa after osmotic stress. *Andrologia.* 2004;36(2):87-93.
124. Fetic S, Yeung CH, Sonntag B, Nieschlag E, Cooper TG. Relationship of cytoplasmic droplets to motility, migration in mucus, and volume regulation of human spermatozoa. *J Androl.* 2006;27(2):294-301.
125. Pelfrey RJ, Overstreet JW, Lewis EL. Abnormalities of sperm morphology in cases of persistent infertility after vasectomy reversal. *Fertil Steril.* 1982;38(1):112-4.
126. Rothmann SA, Bort AM, Quigley J, Pillow R. Sperm morphology classification: a rational method for schemes adopted by the world health organization. *Methods Mol Biol.* 2013;927:27-37.
127. Iwamoto T, Nozawa S, Yoshiike M, Hoshino T, Baba K, Matsushita T et al. Semen quality of 324 fertile Japanese men. *Hum Reprod.* 2006;21(3):760-5.
128. de la Taille A, Rigot JM, Mahe P, Gervais R, Dumur V, Lemaitre L et al. [Correlation of genitourinary abnormalities, spermiogram and CFTR genotype in patients with bilateral agenesis of the vas deferens]. *Prog Urol.* 1998;8(3):370-6.
129. Weiske WH, Salzler N, Schroeder-Printzen I, Weidner W. Clinical findings in congenital absence of the vasa deferentia. *Andrologia.* 2000;32(1):13-8.
130. Daudin M, Bieth E, Bujan L, Massat G, Pontonnier F, Mieusset R. Congenital bilateral absence of the vas deferens: clinical characteristics, biological parameters, cystic fibrosis transmembrane conductance regulator gene mutations, and implications for genetic counseling. *Fertil Steril.* 2000;74(6):1164-74.

131. von Eckardstein S, Cooper TG, Rutscha K, Meschede D, Horst J, Nieschlag E. Seminal plasma characteristics as indicators of cystic fibrosis transmembrane conductance regulator (CFTR) gene mutations in men with obstructive azoospermia. *Fertil Steril*. 2000;73(6):1226-31.
132. Kussler AP, Pimentel AM, Alcoba DD, Liu IP, Brum IS, Capp E et al. Mechanical processing of hyperviscous semen specimens can negatively affect sperm DNA fragmentation. *Int Urol Nephrol*. 2014;46(4):737-42.
133. de la Taille A, Rigot JM, Mahe P, Vankemmel O, Gervais R, Dumur V et al. Correlation between genito-urinary anomalies, semen analysis and CFTR genotype in patients with congenital bilateral absence of the vas deferens. *Br J Urol*. 1998;81(4):614-9.
134. Haugen TB, Grotmol T. pH of human semen. *Int J Androl*. 1998;21(2):105-8.
135. Le Lannou D, Griveau JF, Le Pichon JP, Quero JC. Effects of chamber depth on the motion pattern of human spermatozoa in semen or in capacitating medium. *Hum Reprod*. 1992;7(10):1417-21.
136. Kraemer M, Fillion C, Martin-Pont B, Auger J. Factors influencing human sperm kinematic measurements by the Celltrak computer-assisted sperm analysis system. *Hum Reprod*. 1998;13(3):611-9.
137. Motulsky H. *Intuitive biostatistics*. New York: Oxford University Press; 1995.
138. Seaman EK, Goluboff E, BarChama N, Fisch H. Accuracy of semen counting chambers as determined by the use of latex beads. *Fertil Steril*. 1996;66(4):662-5.
139. Mahmoud AM, Depoorter B, Piens N, Comhaire FH. The performance of 10 different methods for the estimation of sperm concentration. *Fertil Steril*. 1997;68(2):340-5.
140. Björndahl L, Magnusson K, Holmberg R. Validation of Cellvision Disposable 100 Micrometer Hemocytometer With Improved Neubauer Ruling. *Andrology*. 2018;6(S1).
141. Brazil C, Swan SH, Tollner CR, Treece C, Drobnis EZ, Wang C et al. Quality control of laboratory methods for semen evaluation in a multicenter research study. *J Androl*. 2004;25(4):645-56.
142. Douglas-Hamilton DH, Smith NG, Kuster CE, Vermeiden JP, Althouse GC. Capillary-loaded particle fluid dynamics: effect on estimation of sperm concentration. *J Androl*. 2005;26(1):115-22.
143. Douglas-Hamilton DH, Smith NG, Kuster CE, Vermeiden JP, Althouse GC. Particle distribution in low-volume capillary-loaded chambers. *J Androl*. 2005;26(1):107-14.

144. Björndahl L, Barratt CL. Semen analysis: setting standards for the measurement of sperm numbers. *J Androl.* 2005;26(1):11.
145. Jeyendran RS, Van der Ven HH, Perez-Pelaez M, Crabo BG, Zaneveld LJ. Development of an assay to assess the functional integrity of the human sperm membrane and its relationship to other semen characteristics. *J Reprod Fertil.* 1984;70(1):219-28.
146. Hossain AM, Rizk B, Barik S, Huff C, Thorneycroft IH. Time course of hypo-osmotic swellings of human spermatozoa: evidence of ordered transition between swelling subtypes. *Hum Reprod.* 1998;13(6):1578-83.
147. Holmes E. On osmolality and sperm function during processing for assisted reproduction [Academic]. Stockholm, Sweden: Karolinska Institutet; 2020.
148. Meschede D, Keck C, Zander M, Cooper TG, Yeung CH, Nieschlag E. Influence of three different preparation techniques on the results of human sperm morphology analysis. *Int J Androl.* 1993;16(6):362-9.
149. Kruger TF, Ackerman SB, Simmons KF, Swanson RJ, Brugo SS, Acosta AA. A quick, reliable staining technique for human sperm morphology. *Arch Androl.* 1987;18(3):275-7.
150. Coetzee K, Kruger TF, Vandendael A, de Villiers A, Lombard CJ. Comparison of two staining and evaluation methods used for computerized human sperm morphology evaluations. *Andrologia.* 1997;29(3):133-5.
151. Oral E, Yetis O, Elibol F, Senol H, Irez T, Aksu FM. Assessment of human sperm morphology by strict criteria: comparison of wet preparation versus stained with the modified Diff-Quik method. *Arch Androl.* 2002;48(4):307-14.
152. Maree L, du Plessis SS, Menkveld R, van der Horst G. Morphometric dimensions of the human sperm head depend on the staining method used. *Hum Reprod.* 2010;25(6):1369-82.
153. Auger J, Eustache F, Andersen AG, Irvine DS, Jorgensen N, Skakkebaek NE et al. Sperm morphological defects related to environment, lifestyle and medical history of 1001 male partners of pregnant women from four European cities. *Hum Reprod.* 2001;16(12):2710-7.
154. Menkveld RK, T.F. Basic Semen Analysis. In: Acosta AA, editor. Human spermatozoa in assisted reproduction. Baltimore: Williams & Wilkins; 1990. p. 68-84.
155. Aziz N, Buchan I, Taylor C, Kingsland CR, Lewis-Jones I. The sperm deformity index: a reliable predictor of the outcome of oocyte fertilization in vitro. *Fertil Steril.* 1996;66(6):1000-8.

156. Aziz N, Agarwal A, Lewis-Jones I, Sharma RK, Thomas AJ, Jr. Novel associations between specific sperm morphological defects and leukocytospermia. *Fertil Steril*. 2004;82(3):621-7.
157. Aziz N, Saleh RA, Sharma RK, Lewis-Jones I, Esfandiari N, Thomas AJ, Jr. et al. Novel association between sperm reactive oxygen species production, sperm morphological defects, and the sperm deformity index. *Fertil Steril*. 2004;81(2):349-54.
158. David G, Czyglik F, Mayaux MJ, Schwartz D. The success of A.I.D. and semen characteristics: study on 1489 cycles and 192 ejaculates. *Int J Androl*. 1980;3(6):613-9.
159. Auger JE, F. Standardisation de la classification morphologique des spermatozoïdes humains selon la méthode de David modifiée. *Andrologia*. 2000;10:358-73.
160. ESHRE/NAFA. Manual on basic semen analysis (ESHRE Monographs #2). Oxford: Oxford University Press; 2002.
161. Wang WL, Tu CF, Tan YQ. Insight on multiple morphological abnormalities of sperm flagella in male infertility: what is new? *Asian J Androl*. 2020;22(3):236-45.
162. Coutton C, Escoffier J, Martinez G, Arnoult C, Ray PF. Teratozoospermia: spotlight on the main genetic actors in the human. *Hum Reprod Update*. 2015;21(4):455-85.
163. Jørgensen N, Andersen AG, Eustache F, Irvine DS, Suominen J, Petersen JH et al. Regional differences in semen quality in Europe. *Hum Reprod*. 2001;16(5):1012-9.
164. Muratori M, Marchiani S, Tamburrino L, Cambi M, Lotti F, Natali I et al. DNA fragmentation in brighter sperm predicts male fertility independently from age and semen parameters. *Fertil Steril*. 2015;104(3):582-90 e4.
165. Muratori M, Marchiani S, Tamburrino L, Baldi E. Sperm DNA Fragmentation: Mechanisms of Origin. *Adv Exp Med Biol*. 2019;1166:75-85.
166. Sakkas D, Alvarez JG. Sperm DNA fragmentation: mechanisms of origin, impact on reproductive outcome, and analysis. *Fertil Steril*. 2010;93(4):1027-36.
167. Robinson L, Gallos ID, Conner SJ, Rajkhowa M, Miller D, Lewis S et al. The effect of sperm DNA fragmentation on miscarriage rates: a systematic review and meta-analysis. *Hum Reprod*. 2012;27(10):2908-17.
168. Cissen M, Wely MV, Scholten I, Mansell S, Bruin JP, Mol BW et al. Measuring Sperm DNA Fragmentation and Clinical Outcomes of Medically Assisted Reproduction: A Systematic Review and Meta-Analysis. *PLoS One*. 2016;11(11):e0165125.
169. Simon L, Zini A, Dyachenko A, Ciampi A, Carrell DT. A systematic review and meta-analysis to determine the effect of sperm DNA damage on in vitro fertilization and intracytoplasmic sperm injection outcome. *Asian J Androl*. 2017;19(1):80-90.

170. Tan J, Taskin O, Albert A, Bedaiwy MA. Association between sperm DNA fragmentation and idiopathic recurrent pregnancy loss: a systematic review and meta-analysis. *Reprod Biomed Online*. 2019;38(6):951-60.
171. Evgeni E, Charalabopoulos K, Asimakopoulos B. Human sperm DNA fragmentation and its correlation with conventional semen parameters. *J Reprod Infertil*. 2014;15(1):2-14.
172. Robbins DJ, Coleman MS. Initiator role of double stranded DNA in terminal transferase catalyzed polymerization reactions. *Nucleic Acids Res*. 1988;16(7):2943-57.
173. Sharma R, Masaki J, Agarwal A. Sperm DNA fragmentation analysis using the TUNEL assay. *Methods Mol Biol*. 2013;927:121-36.
174. Muratori M, Tamburrino L, Tocci V, Costantino A, Marchiani S, Giachini C et al. Small variations in crucial steps of TUNEL assay coupled to flow cytometry greatly affect measures of sperm DNA fragmentation. *J Androl*. 2010;31(4):336-45.
175. Muratori M, Tamburrino L, Marchiani S, Guido C, Forti G, Baldi E. Critical aspects of detection of sperm DNA fragmentation by TUNEL/flow cytometry. *Syst Biol Reprod Med*. 2010;56(4):277-85.
176. Darzynkiewicz Z, Galkowski D, Zhao H. Analysis of apoptosis by cytometry using TUNEL assay. *Methods*. 2008;44(3):250-4.
177. Dominguez-Fandos D, Camejo MI, Ballesca JL, Oliva R. Human sperm DNA fragmentation: correlation of TUNEL results as assessed by flow cytometry and optical microscopy. *Cytometry A*. 2007;71(12):1011-8.
178. Sergerie M, Laforest G, Bujan L, Bissonnette F, Bleau G. Sperm DNA fragmentation: threshold value in male fertility. *Hum Reprod*. 2005;20(12):3446-51.
179. Sergerie M, Laforest G, Boulanger K, Bissonnette F, Bleau G. Longitudinal study of sperm DNA fragmentation as measured by terminal uridine nick end-labelling assay. *Hum Reprod*. 2005;20(7):1921-7.
180. Antonucci N, Manes S, Corradetti B, Manicardi GC, Borini A, Bizzaro D. A novel in vitro sperm head decondensation protocol for rapid flow cytometric measurement of deoxyribonucleic acid content. *Fertil Steril*. 2013;99(7):1857-61.
181. Mitchell LA, De Iuliis GN, Aitken RJ. The TUNEL assay consistently underestimates DNA damage in human spermatozoa and is influenced by DNA compaction and cell vitality: development of an improved methodology. *Int J Androl*. 2011;34(1):2-13.
182. Ribeiro SC, Muratori M, De Geyter M, De Geyter C. TUNEL labeling with BrdUTP/anti-BrdUTP greatly underestimates the level of sperm DNA fragmentation in semen evaluation. *PLoS One*. 2017;12(8):e0181802.

183. Muratori M, Pellegrino G, Mangone G, Azzari C, Lotti F, Tarozzi N et al. DNA Fragmentation in Viable and Non-Viable Spermatozoa Discriminates Fertile and Subfertile Subjects with Similar Accuracy. *J Clin Med*. 2020;9(5).
184. Martinez MG, Sanchez-Martin P, Dorado-Silva M, Fernandez JL, Girones E, Johnston SD et al. Magnetic-activated cell sorting is not completely effective at reducing sperm DNA fragmentation. *J Assist Reprod Genet*. 2018;35(12):2215-21.
185. Fernandez JL, Cajigal D, Lopez-Fernandez C, Gosalvez J. Assessing sperm DNA fragmentation with the sperm chromatin dispersion test. *Methods Mol Biol*. 2011;682:291-301.
186. Fernandez JL, Muriel L, Rivero MT, Goyanes V, Vazquez R, Alvarez JG. The sperm chromatin dispersion test: a simple method for the determination of sperm DNA fragmentation. *J Androl*. 2003;24(1):59-66.
187. Gosalvez J, Rodriguez-Predreira M, Mosquera A, Lopez-Fernandez C, Esteves SC, Agarwal A et al. Characterisation of a subpopulation of sperm with massive nuclear damage, as recognised with the sperm chromatin dispersion test. *Andrologia*. 2014;46(6):602-9.
188. Klaude M, Eriksson S, Nygren J, Ahnstrom G. The comet assay: mechanisms and technical considerations. *Mutat Res*. 1996;363(2):89-96.
189. McKelvey-Martin VJ, Melia N, Walsh IK, Johnston SR, Hughes CM, Lewis SE et al. Two potential clinical applications of the alkaline single-cell gel electrophoresis assay: (1). Human bladder washings and transitional cell carcinoma of the bladder; and (2). Human sperm and male infertility. *Mutat Res*. 1997;375(2):93-104.
190. Simon L, Carrell DT. Sperm DNA damage measured by comet assay. *Methods Mol Biol*. 2013;927:137-46.
191. Baumgartner A, Cemeli E, Anderson D. The comet assay in male reproductive toxicology. *Cell Biol Toxicol*. 2009;25(1):81-98.
192. Collins AR, Dobson VL, Dusinska M, Kennedy G, Stetina R. The comet assay: what can it really tell us? *Mutat Res*. 1997;375(2):183-93.
193. Konca K, Lankoff A, Banasik A, Lisowska H, Kuszewski T, Gozdz S et al. A cross-platform public domain PC image-analysis program for the comet assay. *Mutat Res*. 2003;534(1-2):15-20.
194. Gonzalez JE, Romero I, Barquinero JF, Garcia O. Automatic analysis of silver-stained comets by CellProfiler software. *Mutat Res*. 2012;748(1-2):60-4.
195. Simon L, Castillo J, Oliva R, Lewis SE. Relationships between human sperm protamines, DNA damage and assisted reproduction outcomes. *Reprod Biomed Online*. 2011;23(6):724-34.

196. Ribas-Maynou J, Garcia-Peiro A, Fernandez-Encinas A, Abad C, Amengual MJ, Prada E et al. Comprehensive analysis of sperm DNA fragmentation by five different assays: TUNEL assay, SCSA, SCD test and alkaline and neutral Comet assay. *Andrology*. 2013;1(5):715-22.
197. Javed A, Talkad MS, Ramaiah MK. Corrigendum to "Evaluation of sperm DNA fragmentation using multiple methods: a comparison of their predictive power for male infertility". *Clin Exp Reprod Med*. 2019;46(4):211.
198. Javed A, Talkad MS, Ramaiah MK. Evaluation of sperm DNA fragmentation using multiple methods: a comparison of their predictive power for male infertility. *Clin Exp Reprod Med*. 2019;46(1):14-21.
199. Evenson DP. Sperm chromatin structure assay (SCSA(R)). *Methods Mol Biol*. 2013;927:147-64.
200. Evenson D, Jost L. Sperm chromatin structure assay is useful for fertility assessment. *Methods Cell Sci*. 2000;22(2-3):169-89.
201. Evenson D, Jost L. Sperm chromatin structure assay for fertility assessment. *Curr Protoc Cytom*. 2001;Chapter 7:Unit 7 13.
202. Evenson DP, Jost LK, Marshall D, Zinaman MJ, Clegg E, Purvis K et al. Utility of the sperm chromatin structure assay as a diagnostic and prognostic tool in the human fertility clinic. *Hum Reprod*. 1999;14(4):1039-49.
203. Giwercman A, Lindstedt L, Larsson M, Bungum M, Spano M, Levine RJ et al. Sperm chromatin structure assay as an independent predictor of fertility in vivo: a case-control study. *Int J Androl*. 2010;33(1):e221-7.
204. Godo A, Blanco J, Vidal F, Sandalinas M, Garcia-Guixé E, Anton E. Altered segregation pattern and numerical chromosome abnormalities interrelate in spermatozoa from Robertsonian translocation carriers. *Reprod Biomed Online*. 2015;31(1):79-88.
205. Wang BN, B.; Tang, D.; Li, R.; Liu, X.; Song, J.; Wang, W.; Liu, Z. . Analysis of Meiotic Segregation Patterns and Interchromosomal Effects in Sperm from 13 Robertsonian Translocations. *Balkan J Med Genet* 2017;20(1):43-50.
206. Enciso M, Alfarawati S, Wells D. Increased numbers of DNA-damaged spermatozoa in samples presenting an elevated rate of numerical chromosome abnormalities. *Hum Reprod*. 2013;28(6):1707-15.
207. Ramasamy R, Scovell JM, Kovac JR, Cook PJ, Lamb DJ, Lipshultz LI. Fluorescence in situ hybridization detects increased sperm aneuploidy in men with recurrent pregnancy loss. *Fertil Steril*. 2015;103(4):906-9 e1.

208. Kohn TP, Kohn JR, Darilek S, Ramasamy R, Lipshultz L. Genetic counseling for men with recurrent pregnancy loss or recurrent implantation failure due to abnormal sperm chromosomal aneuploidy. *J Assist Reprod Genet.* 2016;33(5):571-6.
209. Ramasamy R, Besada S, Lamb DJ. Fluorescent in situ hybridization of human sperm: diagnostics, indications, and therapeutic implications. *Fertil Steril.* 2014;102(6):1534-9.
210. Ryu HM, Lin WW, Lamb DJ, Chuang W, Lipshultz LI, Bischoff FZ. Increased chromosome X, Y, and 18 nondisjunction in sperm from infertile patients that were identified as normal by strict morphology: implication for intracytoplasmic sperm injection. *Fertil Steril.* 2001;76(5):879-83.
211. Ioannou D, Fortun J, Tempest HG. Meiotic nondisjunction and sperm aneuploidy in humans. *Reproduction.* 2019;157(1):R15-R31.
212. Neusser M, Rogenhofer N, Durl S, Ochsenkuhn R, Trottman M, Jurinovic V et al. Increased chromosome 16 disomy rates in human spermatozoa and recurrent spontaneous abortions. *Fertil Steril.* 2015;104(5):1130-7 e1-10.
213. Templado C, Uroz L, Estop A. New insights on the origin and relevance of aneuploidy in human spermatozoa. *Mol Hum Reprod.* 2013;19(10):634-43.
214. Tomlinson MJ, Barratt CL, Cooke ID. Prospective study of leukocytes and leukocyte subpopulations in semen suggests they are not a cause of male infertility. *Fertil Steril.* 1993;60(6):1069-75.
215. Wolff H. The biologic significance of white blood cells in semen. *Fertil Steril.* 1995;63(6):1143-57.
216. Homyk M, Anderson DJ, Wolff H, Herr JC. Differential diagnosis of immature germ cells in semen utilizing monoclonal antibody MHS-10 to the intra-acrosomal antigen SP-10. *Fertil Steril.* 1990;53(2):323-30.
217. Nahoum CR, Cardozo D. Staining for volumetric count of leukocytes in semen and prostrate-vesicular fluid. *Fertil Steril.* 1980;34(1):68-9.
218. Sharma RK, Pasqualotto AE, Nelson DR, Thomas AJ, Jr., Agarwal A. Relationship between seminal white blood cell counts and oxidative stress in men treated at an infertility clinic. *J Androl.* 2001;22(4):575-83.
219. Punab M, Loivukene K, Kermes K, Mandar R. The limit of leucocytospermia from the microbiological viewpoint. *Andrologia.* 2003;35(5):271-8.
220. Aitken RJ, Baker HW. Seminal leukocytes: passengers, terrorists or good samaritans? *Hum Reprod.* 1995;10(7):1736-9.

221. Rossi AG, Aitken RJ. Interactions between leukocytes and the male reproductive system. The unanswered questions. *Adv Exp Med Biol.* 1997;424:245-52.
222. Couture M, Ulstein M, Leonard J, Paulsen CA. Improved staining method for differentiating immature germ cells from white blood cells in human seminal fluid. *Andrologia.* 1976;8(1):61-6.
223. Ezeh UI, Martin M, Cooke ID, Moore HD. Correlation of testicular pathology and sperm extraction in azoospermic men with ejaculated spermatids detected by immunofluorescent localization. *Hum Reprod.* 1998;13(11):3061-5.
224. Barbonetti A, Castellini C, D'Andrea S, Cordeschi G, Santucci R, Francavilla S et al. Prevalence of anti-sperm antibodies and relationship of degree of sperm auto-immunization to semen parameters and post-coital test outcome: a retrospective analysis of over 10 000 men. *Hum Reprod.* 2019;34(5):834-41.
225. Kremer J, Jager S. Characteristics of anti-spermatozoal antibodies responsible for the shaking phenomenon with special regard to immunoglobulin class and antigen-reactive sites. *Int J Androl.* 1980;3(2):143-52.
226. Bronson R, Cooper G, Rosenfeld D. Sperm antibodies: their role in infertility. *Fertil Steril.* 1984;42(2):171-83.
227. Bronson R, Cooper G, Rosenfeld D, Witkin SS. Detection of spontaneously occurring sperm-directed antibodies in infertile couples by immunobead binding and enzyme-linked immunosorbent assay. *Ann N Y Acad Sci.* 1984;438:504-7.
228. Bronson RA. Detection of sperm specific antibodies on the spermatozoa surface by immunobead binding. *Arch Androl.* 1982;9:61.
229. Clarke GN. Immunoglobulin class of sperm-bound antibodies in semen. . In: Bratanov K, editor. *Immunology of Reproduction.* Blgarian: Bulgarian Academy of Sciences Press; 1982. p. 462-85.
230. Hellstrom WJ, Samuels SJ, Waits AB, Overstreet JW. A comparison of the usefulness of SpermMar and immunobead tests for the detection of antisperm antibodies. *Fertil Steril.* 1989;52(6):1027-31.
231. Scarselli G, Livi C, Chelo E, Dubini V, Pellegrini S. Approach to immunological male infertility: a comparison between MAR test and direct immunobead test. *Acta Eur Fertil.* 1987;18(1):55-7.
232. MacMillan RA, Baker HW. Comparison of latex and polyacrylamide beads for detecting sperm antibodies. *Clin Reprod Fertil.* 1987;5(4):203-9.
233. Meinertz H, Bronson R. Detection of antisperm antibodies on the surface of motile spermatozoa. Comparison of the immunobead binding technique (IBT) and the mixed antiglobulin reaction (MAR). *Am J Reprod Immunol Microbiol.* 1988;18(4):120-3.

234. Abshagen K, Behre HM, Cooper TG, Nieschlag E. Influence of sperm surface antibodies on spontaneous pregnancy rates. *Fertil Steril*. 1998;70(2):355-6.
235. Chiu WW, Chamley LW. Clinical associations and mechanisms of action of antisperm antibodies. *Fertil Steril*. 2004;82(3):529-35.
236. Rajah SV, Parslow JM, Howell RJ, Hendry WF. Comparison of mixed antiglobulin reaction and direct immunobead test for detection of sperm-bound antibodies in subfertile males. *Fertil Steril*. 1992;57(6):1300-3.
237. Gould JE, Brazil CK, Overstreet JW. Sperm-immunobead binding decreases with in vitro incubation. *Fertil Steril*. 1994;62(1):167-71.
238. Barratt CL, Dunphy BC, McLeod I, Cooke ID. The poor prognostic value of low to moderate levels of sperm surface-bound antibodies. *Hum Reprod*. 1992;7(1):95-8.
239. von der Kammer H, Scheit KH, Weidner W, Cooper TG. The evaluation of markers of prostatic function. *Urol Res*. 1991;19(6):343-7.
240. Cooper TG, Weidner W, Nieschlag E. The influence of inflammation of the human male genital tract on secretion of the seminal markers alpha-glucosidase, glycerophosphocholine, carnitine, fructose and citric acid. *Int J Androl*. 1990;13(5):329-36.
241. Moellering H, Gruber W. Determination of citrate with citrate lyase. *Anal Biochem*. 1966;17(3):369-76.
242. Heite HJ, Wetterauer W. [Acid phosphatase in seminal fluid--method of estimation and diagnostic significance]. *Andrologia*. 1979;11(2):113-22.
243. Wang J, Niu Y, Zhang C, Chen Y. A micro-plate colorimetric assay for rapid determination of trace zinc in animal feed, pet food and drinking water by ion masking and statistical partitioning correction. *Food Chem*. 2018;245:337-45.
244. Johnsen O, Eliasson R. Evaluation of a commercially available kit for the colorimetric determination of zinc in human seminal plasma. *Int J Androl*. 1987;10(2):435-40.
245. Cooper TG, Jockenhovel F, Nieschlag E. Variations in semen parameters from fathers. *Hum Reprod*. 1991;6(6):859-66.
246. Björndahl L. Prevalence of high zinc concentrations in 45,000 ejaculates - Unpublished data. 2021.
247. Karvonen MJ, Malm M. Colorimetric determination of fructose with indol. *Scand J Clin Lab Invest*. 1955;7(4):305-7.

248. Paquin R, Chapdelaine P, Dube JY, Tremblay RR. Similar biochemical properties of human seminal plasma and epididymal alpha-1,4-glucosidase. *J Androl.* 1984;5(4):277-82.
249. Cooper TG, Yeung CH, Nashan D, Jockenhovel F, Nieschlag E. Improvement in the assessment of human epididymal function by the use of inhibitors in the assay of alpha-glucosidase in seminal plasma. *Int J Androl.* 1990;13(4):297-305.
250. Björndahl L, Magnusson K, Larsson Chatziantonis M, Holmberg R, Flanagan J. When is a vasectomy successful? Laboratory aspects. *Andrology.* 2019;7, Suppl 1:96-97.
251. MacLeod J, Hotchkiss RS. The Distribution of Spermatozoa and of Certain Chemical Constituents in the Human Ejaculate. *J Urol.* 1942;48(2):225-9.
252. MacLeod J, Gold RZ. The male factor in fertility and infertility. III. An analysis of motile activity in the spermatozoa of 1000 fertile men and 1000 men in infertile marriage. *Fertil Steril.* 1951;2(3):187-204.
253. Björndahl L, Kvist U. Influence of seminal vesicular fluid on the zinc content of human sperm chromatin. *Int J Androl.* 1990;13(3):232-7.
254. Kvist U, Kjellberg S, Björndahl L, Soufir JC, Arver S. Seminal fluid from men with agenesis of the Wolffian ducts: zinc-binding properties and effects on sperm chromatin stability. *Int J Androl.* 1990;13(4):245-52.
255. Björndahl L, Kjellberg S, Kvist U. Ejaculatory sequence in men with low sperm chromatin-zinc. *Int J Androl.* 1991;14(3):174-8.
256. Lindholmer C. Survival of human spermatozoa in different fractions of split ejaculate. *Fertil Steril.* 1973;24(7):521-6.
257. Lundquist F. Aspects of the biochemistry of human semen. . *Acta Physiol Scand.* 1949;19, Suppl 66:1-108.
258. Beiswanger JC, Deaton JL, Jarow JP. Partial ejaculatory duct obstruction causing early demise of sperm. *Urology.* 1998;51(1):125-7.
259. Eisenberg ML, Walsh TJ, Garcia MM, Shinohara K, Turek PJ. Ejaculatory duct manometry in normal men and in patients with ejaculatory duct obstruction. *J Urol.* 2008;180(1):255-60; discussion 60.
260. Fisch H, Kang YM, Johnson CW, Goluboff ET. Ejaculatory duct obstruction. *Curr Opin Urol.* 2002;12(6):509-15.
261. Lotti F, Corona G, Cocci A, Cipriani S, Baldi E, Degl'Innocenti S et al. The prevalence of midline prostatic cysts and the relationship between cyst size and semen parameters among infertile and fertile men. *Hum Reprod.* 2018;33(11):2023-34.

262. Nagler HM, Rotman M, Zoltan E, Fisch H. The natural history of partial ejaculatory duct obstruction. *J Urol*. 2002;167(1):253-4.
263. Onur MR, Orhan I, Firdolas F, Onur R, Kocakoc E. Clinical and radiological evaluation of ejaculatory duct obstruction. *Arch Androl*. 2007;53(4):179-86.
264. Smith JF, Walsh TJ, Turek PJ. Ejaculatory duct obstruction. *Urol Clin North Am*. 2008;35(2):221-7, viii.
265. Turek PJ, Magana JO, Lipshultz LI. Semen parameters before and after transurethral surgery for ejaculatory duct obstruction. *J Urol*. 1996;155(4):1291-3.
266. Avellino GJ, Lipshultz LI, Sigman M, Hwang K. Transurethral resection of the ejaculatory ducts: etiology of obstruction and surgical treatment options. *Fertil Steril*. 2019;111(3):427-43.
267. Chan PTK, Aube M. Infertility due to ejaculatory duct obstruction can be surgically corrected-a forgotten lesson in male reproductive urology. *Fertil Steril*. 2018;110(7):1275-6.
268. El-Assmy A, El-Tholoth H, Abouelkheir RT, Abou-El-Ghar ME. Transurethral resection of ejaculatory duct in infertile men: outcome and predictors of success. *Int Urol Nephrol*. 2012;44(6):1623-30.
269. Halpern JA, Brannigan RE, Schlegel PN. Fertility-enhancing male reproductive surgery: glimpses into the past and thoughts for the future. *Fertil Steril*. 2019;112(3):426-37.
270. Kadioglu A, Cayan S, Tefekli A, Orhan I, Engin G, Turek PJ. Does response to treatment of ejaculatory duct obstruction in infertile men vary with pathology? *Fertil Steril*. 2001;76(1):138-42.
271. McQuaid JW, Tanrikut C. Ejaculatory duct obstruction: current diagnosis and treatment. *Curr Urol Rep*. 2013;14(4):291-7.
272. Fisch H, Lambert SM, Goluboff ET. Management of ejaculatory duct obstruction: etiology, diagnosis, and treatment. *World J Urol*. 2006;24(6):604-10.
273. Mortimer D, Björndahl L, Barratt CLR, Castilla JA, Menkveld R, Kvist U et al. A practical guide to basic laboratory andrology. 1st ed. Cambridge, UK ; New York: Cambridge University Press; In Press.
274. Amelar RD, Hotchkiss RS. The Split Ejaculate: Its Use in the Management of Male Infertility. *Fertil Steril*. 1965;16:46-60.
275. Eisenberg ML, Li S, Behr B, Pera RR, Cullen MR. Relationship between semen production and medical comorbidity. *Fertil Steril*. 2015;103(1):66-71.

276. Selvaraju V, Baskaran S, Agarwal A, Henkel R. Environmental contaminants and male infertility: Effects and mechanisms. *Andrologia*. 2020:e13646.
277. Levine H, Jorgensen N, Martino-Andrade A, Mendiola J, Weksler-Derri D, Mindlis I et al. Temporal trends in sperm count: a systematic review and meta-regression analysis. *Hum Reprod Update*. 2017;23(6):646-59.
278. Muratori M, De Geyter C. Chromatin condensation, fragmentation of DNA and differences in the epigenetic signature of infertile men. *Best Pract Res Clin Endocrinol Metab*. 2019;33(1):117-26.
279. Jones R, Mann T, Sherins R. Peroxidative breakdown of phospholipids in human spermatozoa, spermicidal properties of fatty acid peroxides, and protective action of seminal plasma. *Fertil Steril*. 1979;31(5):531-7.
280. Alvarez JG, Touchstone JC, Blasco L, Storey BT. Spontaneous lipid peroxidation and production of hydrogen peroxide and superoxide in human spermatozoa. Superoxide dismutase as major enzyme protectant against oxygen toxicity. *J Androl*. 1987;8(5):338-48.
281. Aitken RJ, Clarkson JS. Cellular basis of defective sperm function and its association with the genesis of reactive oxygen species by human spermatozoa. *J Reprod Fertil*. 1987;81(2):459-69.
282. Bisht S, Dada R. Oxidative stress: Major executioner in disease pathology, role in sperm DNA damage and preventive strategies. *Front Biosci (Schol Ed)*. 2017;9:420-47.
283. O'Flaherty C, Matsushita-Fournier D. Reactive oxygen species and protein modifications in spermatozoa. *Biol Reprod*. 2017;97(4):577-85.
284. Aitken RJ. Impact of oxidative stress on male and female germ cells: implications for fertility. *Reproduction*. 2020;159(4):R189-R201.
285. Dias TR. Measurement of Reactive Oxygen Species in Semen Samples Using Chemiluminescence. *Methods Mol Biol*. 2021;2202:103-9.
286. Miller NJ, Rice-Evans C, Davies MJ. A new method for measuring antioxidant activity. *Biochem Soc Trans*. 1993;21(2):95S.
287. Miller NJ, Rice-Evans CA. Factors influencing the antioxidant activity determined by the ABTS.+ radical cation assay. *Free Radic Res*. 1997;26(3):195-9.
288. Xu F, Guo G, Zhu W, Fan L. Human sperm acrosome function assays are predictive of fertilization rate in vitro: a retrospective cohort study and meta-analysis. *Reprod Biol Endocrinol*. 2018;16(1):81.

289. Franken DR, Bastiaan HS, Oehninger SC. Physiological induction of the acrosome reaction in human sperm: validation of a microassay using minimal volumes of solubilized, homologous zona pellucida. *J Assist Reprod Genet.* 2000;17(3):156-61.
290. Tamburrino L, Marchiani S, Muratori M, Luconi M, Baldi E. Progesterone, spermatozoa and reproduction: An updated review. *Mol Cell Endocrinol.* 2020;516:110952.
291. Liu DY, Baker HW. Relationship between the zona pellucida (ZP) and ionophore A23187-induced acrosome reaction and the ability of sperm to penetrate the ZP in men with normal sperm-ZP binding. *Fertil Steril.* 1996;66(2):312-5.
292. Fierro R, Foliguet B, Grignon G, Daniel M, Bene MC, Faure GC et al. Lectin-binding sites on human sperm during acrosome reaction: modifications judged by electron microscopy/flow cytometry. *Arch Androl.* 1996;36(3):187-96.
293. Cross N. Methods for evaluating the acrosomal status of human sperm. In: Fenichel P, Parinaud J, editors. Human sperm acrosome reaction. Colloques INSERM. Paris: John Libbey Eurotext; 1995. p. 277-85.
294. Aitken RJ, Buckingham DW, Fang HG. Analysis of the responses of human spermatozoa to A23187 employing a novel technique for assessing the acrosome reaction. *J Androl.* 1993;14(2):132-41.
295. Jin SK, Yang WX. Factors and pathways involved in capacitation: how are they regulated? *Oncotarget.* 2017;8(2):3600-27.
296. Cross N. Two simple methods for detecting acrosome-reacted human sperm. *Gamete Research.* 1986;15 213-26.
297. Mortimer D, Curtis EF, Miller RG. Specific labelling by peanut agglutinin of the outer acrosomal membrane of the human spermatozoon. *J Reprod Fertil.* 1987;81(1):127-35.
298. Pember S, Jequier A, Yovich J. A Test of the Human Sperm Acrosome Reaction Following Ionophore Challenge Relationship to Fertility and Other Seminal Parameters. *Journal of andrology.* 1991;12:98-103.
299. Krausz C, Bonaccorsi L, Maggio P, Luconi M, Criscuoli L, Fuzzi B et al. Two functional assays of sperm responsiveness to progesterone and their predictive values in in-vitro fertilization. *Hum Reprod.* 1996;11(8):1661-7.
300. Björndahl L, Kvist U. A model for the importance of zinc in the dynamics of human sperm chromatin stabilization after ejaculation in relation to sperm DNA vulnerability. *Syst Biol Reprod Med.* 2011;57(1-2):86-92.
301. Colaco S, Sakkas D. Paternal factors contributing to embryo quality. *J Assist Reprod Genet.* 2018;35(11):1953-68.

302. Marchiani S, Tamburrino L, Benini F, Fanfani L, Dolce R, Rastrelli G et al. Chromatin Protamination and Catsper Expression in Spermatozoa Predict Clinical Outcomes after Assisted Reproduction Programs. *Sci Rep*. 2017;7(1):15122.
303. Marchiani S, Tamburrino L, Olivito B, Betti L, Azzari C, Forti G et al. Characterization and sorting of flow cytometric populations in human semen. *Andrology*. 2014;2(3):394-401.
304. Smith JF, Syritsyna O, Fellous M, Serres C, Mannowetz N, Kirichok Y et al. Disruption of the principal, progesterone-activated sperm Ca^{2+} channel in a CatSper2-deficient infertile patient. *Proc Natl Acad Sci U S A*. 2013;110(17):6823-8.
305. Ren D, Navarro B, Perez G, Jackson AC, Hsu S, Shi Q et al. A sperm ion channel required for sperm motility and male fertility. *Nature*. 2001;413(6856):603-9.
306. Brown SG, Publicover SJ, Barratt CLR, Martins da Silva SJ. Human sperm ion channel (dys)function: implications for fertilization. *Hum Reprod Update*. 2019;25(6):758-76.
307. Luo T, Chen HY, Zou QX, Wang T, Cheng YM, Wang HF et al. A novel copy number variation in CATSPER2 causes idiopathic male infertility with normal semen parameters. *Hum Reprod*. 2019;34(3):414-23.
308. Kelly MC, Brown SG, Costello SM, Ramalingam M, Drew E, Publicover SJ et al. Single-cell analysis of $[\text{Ca}^{2+}]_i$ signalling in sub-fertile men: characteristics and relation to fertilization outcome. *Hum Reprod*. 2018;33(6):1023-33.
309. Davis RO, Katz DF. Standardization and comparability of CASA instruments. *J Androl*. 1992;13(1):81-6.
310. Mortimer D, Mortimer S, editors. Routine application of CASA in human clinical andrology and ART laboratories. 13th International Symposium on Spermatology; 2018; Stockholm New York: Springer; 2021.
311. Mortimer D, Aitken RJ, Mortimer ST, Pacey AA. Workshop report: clinical CASA--the quest for consensus. *Reproduction, fertility, and development*. 1995;7(4):951-9.
312. ESHRE Special Interest Group in Andrology. Guidelines on the application of CASA technology in the analysis of spermatozoa. 1998 Jan. Report No.: 0268-1161 (Print) Contract No.: 13.
313. Freour T, Jean M, Mirallie S, Barriere P. Computer-assisted sperm analysis parameters in young fertile sperm donors and relationship with age. *Syst Biol Reprod Med*. 2012;58(2):102-6.
314. Yuan P, Zhang Z, Luo C, Quan S, Chu Q, Zhao S. [Computer-assisted sperm analysis for assessing sperm mobility parameters in in vitro fertilization]. *Nan Fang Yi Ke Da Xue Xue Bao*. 2013;33(3):448-50.

315. Cupples G, Gallagher MT, Smith AU, Kirkman-Brown J, editors. Heads and tails. Requirements for informative and robust computational measures of sperm motility. 13th International Symposium on Spermatology; 2018; Stockholm. New York: Springer; 2021.
316. Gallagher MT, Cupples G, Ooi EH, Kirkman-Brown JC, Smith DJ. Rapid sperm capture: high-throughput flagellar waveform analysis. *Hum Reprod.* 2019;34(7):1173-85.
317. Mortimer ST. A critical review of the physiological importance and analysis of sperm movement in mammals. *Hum Reprod Update.* 1997;3(5):403-39.
318. Zhu JJ, Pacey AA, Barratt CL, Cooke ID. Computer-assisted measurement of hyperactivation in human spermatozoa: differences between European and American versions of the Hamilton-Thorn motility analyser. *Hum Reprod.* 1994;9(3):456-62.
319. Mortimer D, Mortimer ST. Computer-Aided Sperm Analysis (CASA) of sperm motility and hyperactivation. *Methods Mol Biol.* 2013;927:77-87.
320. Goodson SG, White S, Stevans AM, Bhat S, Kao CY, Jaworski S et al. CASAnova: a multiclass support vector machine model for the classification of human sperm motility patterns. *Biol Reprod.* 2017;97(5):698-708.
321. Ooi EH, Smith DJ, Gadelha H, Gaffney EA, Kirkman-Brown J. The mechanics of hyperactivation in adhered human sperm. *R Soc Open Sci.* 2014;1(2):140230.
322. Garrett C, Baker HW. A new fully automated system for the morphometric analysis of human sperm heads. *Fertil Steril.* 1995;63(6):1306-17.
323. Lacquet FA, Kruger TF, Du Toit TC, Lombard CJ, Sanchez Sarmiento CA, De Villiers A et al. Slide preparation and staining procedures for reliable results using computerized morphology. *Arch Androl.* 1996;36(2):133-8.
324. Menkveld R, Lacquet FA, Kruger TF, Lombard CJ, Sanchez Sarmiento CA, de Villiers A. Effects of different staining and washing procedures on the results of human sperm morphology evaluation by manual and computerised methods. *Andrologia.* 1997;29(1):1-7.
325. Coetzee K, Kruger TF, Lombard CJ. Repeatability and variance analysis on multiple computer-assisted (IVOS) sperm morphology readings. *Andrologia.* 1999;31(3):163-8.
326. Coetzee K, Kruger TF, Lombard CJ, Shaughnessy D, Oehninger S, Ozgur K et al. Assessment of interlaboratory and intralaboratory sperm morphology readings with the use of a Hamilton Thorne Research integrated visual optical system semen analyzer. *Fertil Steril.* 1999;71(1):80-4.

327. Coetzee K, de Villiers A, Kruger TF, Lombard CJ. Clinical value of using an automated sperm morphology analyzer (IVOS). *Fertil Steril*. 1999;71(2):222-5.
328. An LM, Ma JH, Li WW, Xia XY, Huang YF, Wang WB. [Flow cytometry for evaluation of sperm plasma membrane integrity in varicocele patients]. *Zhonghua Nan Ke Xue*. 2011;17(11):984-8.
329. Isobe T. New method to estimate the possibility of natural pregnancy using computer-assisted sperm analysis. *Syst Biol Reprod Med*. 2012;58(6):339-47.
330. Yaniz JL, Vicente-Fiel S, Soler C, Recreo P, Carretero T, Bono A et al. Comparison of different statistical approaches to evaluate morphometric sperm subpopulations in men. *Asian J Androl*. 2016;18(6):819-23.
331. Lu JC, Huang YF, Lu NQ. Computer-aided sperm analysis: past, present and future. *Andrologia*. 2014;46(4):329-38.
332. Vested A, Ramlau-Hansen CH, Bonde JP, Thulstrup AM, Kristensen SL, Toft G. A comparison of conventional and computer-assisted semen analysis (CRISMAS software) using samples from 166 young Danish men. *Asian J Androl*. 2011;13(3):453-8.
333. Hu YA, Lu JC, Shao Y, Huang YF, Lu NQ. Comparison of the semen analysis results obtained from two branded computer-aided sperm analysis systems. *Andrologia*. 2013;45(5):315-8.
334. Dearing C, Jayasena C, Lindsay K. Can the Sperm Class Analyser (SCA) CASA-Mot system for human sperm motility analysis reduce imprecision and operator subjectivity and improve semen analysis? *Hum Fertil (Camb)*. 2019:1-11.
335. Dearing CG, Kilburn S, Lindsay KS. Validation of the sperm class analyser CASA system for sperm counting in a busy diagnostic semen analysis laboratory. *Hum Fertil (Camb)*. 2014;17(1):37-44.
336. Kose M, Sokmensuer LK, Demir A, Bozdog G, Gunalp S. Manual versus computer-automated semen analysis. *Clin Exp Obstet Gynecol*. 2014;41(6):662-4.
337. Lammers J, Splingart C, Barriere P, Jean M, Freour T. Double-blind prospective study comparing two automated sperm analyzers versus manual semen assessment. *J Assist Reprod Genet*. 2014;31(1):35-43.
338. Urbano LF, Masson P, VerMilyea M, Kam M. Automatic Tracking and Motility Analysis of Human Sperm in Time-Lapse Images. *IEEE Trans Med Imaging*. 2017;36(3):792-801.
339. Wei SY, Chao HH, Huang HP, Hsu CF, Li SH, Hsu L. A collective tracking method for preliminary sperm analysis. *Biomed Eng Online*. 2019;18(1):112.

340. Elsayed M, El-Sherry TM, Abdelgawad M. Development of computer-assisted sperm analysis plugin for analyzing sperm motion in microfluidic environments using Image-J. *Theriogenology*. 2015;84(8):1367-77.
341. Kobori Y, Pfanner P, Prins GS, Niederberger C. Novel device for male infertility screening with single-ball lens microscope and smartphone. *Fertil Steril*. 2016;106(3):574-8.
342. Holmes E, Bjorndahl L, Kvist U. Hypotonic challenge reduces human sperm motility through coiling and folding of the tail. *Andrologia*. 2020:e13859.
343. Mortimer ST, Swan MA, Mortimer D. Effect of seminal plasma on capacitation and hyperactivation in human spermatozoa. *Hum Reprod*. 1998;13(8):2139-46.
344. Canale D, Giorgi PM, Gasperini M, Pucci E, Barletta D, Gasperi M et al. Inter and intra-individual variability of sperm morphology after selection with three different techniques: layering, swimup from pellet and percoll. *J Endocrinol Invest*. 1994;17(9):729-32.
345. Henkel R. Sperm preparation: state-of-the-art--physiological aspects and application of advanced sperm preparation methods. *Asian J Androl*. 2012;14(2):260-9.
346. Volpes A, Sammartano F, Rizzari S, Gullo S, Marino A, Allegra A. The pellet swim-up is the best technique for sperm preparation during in vitro fertilization procedures. *J Assist Reprod Genet*. 2016;33(6):765-70.
347. Henkel RR, Schill WB. Sperm preparation for ART. *Reprod Biol Endocrinol*. 2003;1:108.
348. Aitken RJ, Clarkson JS. Significance of reactive oxygen species and antioxidants in defining the efficacy of sperm preparation techniques. *J Androl*. 1988;9(6):367-76.
349. Iwasaki A, Gagnon C. Formation of reactive oxygen species in spermatozoa of infertile patients. *Fertil Steril*. 1992;57(2):409-16.
350. Boomsma CM, Heineman MJ, Cohlen BJ, Farquhar C. Semen preparation techniques for intrauterine insemination. *Cochrane Database Syst Rev*. 2007(4):Cd004507.
351. Simon L, Emery B, Carrell DT. Sperm DNA Fragmentation: Consequences for Reproduction. *Adv Exp Med Biol*. 2019;1166:87-105.
352. Jayaraman V, Upadhy D, Narayan PK, Adiga SK. Sperm processing by swim-up and density gradient is effective in elimination of sperm with DNA damage. *J Assist Reprod Genet*. 2012;29(6):557-63.

353. Aitken RJ, Finnie JM, Muscio L, Whiting S, Connaughton HS, Kuczera L et al. Potential importance of transition metals in the induction of DNA damage by sperm preparation media. *Hum Reprod.* 2014;29(10):2136-47.
354. Muratori M, Tarozzi N, Cambi M, Boni L, Iorio AL, Passaro C et al. Variation of DNA Fragmentation Levels During Density Gradient Sperm Selection for Assisted Reproduction Techniques: A Possible New Male Predictive Parameter of Pregnancy? *Medicine (Baltimore).* 2016;95(20):e3624.
355. Björndahl L, Mohammadieh M, Pourian M, Soderlund I, Kvist U. Contamination by seminal plasma factors during sperm selection. *J Androl.* 2005;26(2):170-3.
356. Oseguera-Lopez I, Ruiz-Diaz S, Ramos-Ibeas P, Perez-Cerezales S. Novel Techniques of Sperm Selection for Improving IVF and ICSI Outcomes. *Front Cell Dev Biol.* 2019;7:298.
357. International Standards Organization. ISO 15189:2012 Medical Laboratories – Requirements for Quality and Competence. Geneva 2012.
358. Lepine S, McDowell S, Searle LM, Kroon B, Glujovsky D, Yazdani A. Advanced sperm selection techniques for assisted reproduction. *Cochrane Database Syst Rev.* 2019;7:CD010461.
359. Bucar S, Goncalves A, Rocha E, Barros A, Sousa M, Sa R. DNA fragmentation in human sperm after magnetic-activated cell sorting. *J Assist Reprod Genet.* 2015;32(1):147-54.
360. Zafer M, Horvath H, Mmeje O, van der Poel S, Semprini AE, Rutherford G et al. Effectiveness of semen washing to prevent human immunodeficiency virus (HIV) transmission and assist pregnancy in HIV-discordant couples: a systematic review and meta-analysis. *Fertil Steril.* 2016;105(3):645-55 e2.
361. Gilling-Smith C, Nicopoullos JD, Semprini AE, Frodsham LC. HIV and reproductive care--a review of current practice. *BJOG.* 2006;113(8):869-78.
362. Savasi V, Ferrazzi E, Lanzani C, Oneta M, Parrilla B, Persico T. Safety of sperm washing and ART outcome in 741 HIV-1-serodiscordant couples. *Hum Reprod.* 2007;22(3):772-7.
363. Fourie JM, Loskutoff N, Huyser C. Semen decontamination for the elimination of seminal HIV-1. *Reprod Biomed Online.* 2015;30(3):296-302.
364. Gilling-Smith C, Emiliani S, Almeida P, Liesnard C, Englert Y. Laboratory safety during assisted reproduction in patients with blood-borne viruses. *Hum Reprod.* 2005;20(6):1433-8.

365. Loutfy MR, Wu W, Letchumanan M, Bondy L, Antoniou T, Margolese S et al. Systematic review of HIV transmission between heterosexual serodiscordant couples where the HIV-positive partner is fully suppressed on antiretroviral therapy. *PLoS One*. 2013;8(2):e55747.
366. LeMessurier J, Traversy G, Varsaneux O, Weekes M, Avey MT, Niragira O et al. Risk of sexual transmission of human immunodeficiency virus with antiretroviral therapy, suppressed viral load and condom use: a systematic review. *CMAJ*. 2018;190(46):E1350-E60.
367. Rodger AJ, Cambiano V, Bruun T, Vernazza P, Collins S, van Lunzen J et al. Sexual Activity Without Condoms and Risk of HIV Transmission in Serodifferent Couples When the HIV-Positive Partner Is Using Suppressive Antiretroviral Therapy. *JAMA*. 2016;316(2):171-81.
368. Del Romero J, Rio I, Castilla J, Baza B, Paredes V, Vera M et al. Absence of transmission from HIV-infected individuals with HAART to their heterosexual serodiscordant partners. *Enferm Infecc Microbiol Clin*. 2015;33(10):666-72.
369. Mahadevan M, Leeton JF, Trounson AO. Noninvasive method of semen collection for successful artificial insemination in a case of retrograde ejaculation. *Fertil Steril*. 1981;36(2):243-7.
370. Bunge RG, Keettel WC, Sherman JK. Clinical use of frozen semen: report of four cases. *Fertil Steril*. 1954;5(6):520-9.
371. Bunge RG, Sherman JK. Fertilizing capacity of frozen human spermatozoa. *Nature*. 1953;172(4382):767-8.
372. Polge C, Smith AU, Parkes AS. Revival of spermatozoa after vitrification and dehydration at low temperatures. *Nature*. 1949;164(4172):666.
373. Clarke GN, Bourne H, Hill P, Johnston WI, Speirs A, McBain JC et al. Artificial insemination and in-vitro fertilization using donor spermatozoa: a report on 15 years of experience. *Hum Reprod*. 1997;12(4):722-6.
374. David G, Czyglik F, Mayaux MJ, Martin-Boyce A, Schwartz D. Artificial insemination with frozen sperm: protocol, method of analysis and results for 1188 women. *Br J Obstet Gynaecol*. 1980;87(11):1022-8.
375. Leibo SP, Songsasen N. Cryopreservation of gametes and embryos of non-domestic species. *Theriogenology*. 2002;57(1):303-26.
376. Perloff WH, Steinberger E, Sherman JK. Conception with Human Spermatozoa Frozen by Nitrogen Vapor Technic. *Fertil Steril*. 1964;15:501-4.
377. Mortimer D. Current and future concepts and practices in human sperm cryobanking. *Reprod Biomed Online*. 2004;9(2):134-51.

378. Keel BA, Webster BW. Semen cryopreservation methodology and results. In: Barratt CLR, Cooke ID, editors. Donor insemination. Cambridge: Cambridge University Press; 1993. p. 71-96.
379. Sherman JK. Cryopreservation of human semen. In: Keel BA, Webster BW, editors. CRC handbook of the laboratory diagnosis and treatment of infertility. Boca Raton: CRC Press 1990. p. 229-59.
380. Watson PF. Recent developments and concepts in the cryopreservation of spermatozoa and the assessment of their post-thawing function. *Reproduction, fertility, and development*. 1995;7(4):871-91.
381. Clarke GN, Liu DY, Baker HW. Improved sperm cryopreservation using cold cryoprotectant. *Reproduction, fertility, and development*. 2003;15(7-8):377-81.
382. Nijs M, Creemers E, Cox A, Janssen M, Vanheusden E, Castro-Sanchez Y et al. Influence of freeze-thawing on hyaluronic acid binding of human spermatozoa. *Reprod Biomed Online*. 2009;19(2):202-6.
383. Guan HT, Zheng Y, Wang JJ, Meng TQ, Xia W, Hu SH et al. Relationship between donor sperm parameters and pregnancy outcome after intrauterine insemination: analysis of 2821 cycles in 1355 couples. *Andrologia*. 2016;48(1):29-36.
384. Clarke GN, Liu DY, Baker HW. Recovery of human sperm motility and ability to interact with the human zona pellucida after more than 28 years of storage in liquid nitrogen. *Fertil Steril*. 2006;86(3):721-2.
385. Feldschuh J, Brassel J, Durso N, Levine A. Successful sperm storage for 28 years. *Fertil Steril*. 2005;84(4):1017.
386. Sharma R, Kattoor AJ, Ghulmiyyah J, Agarwal A. Effect of sperm storage and selection techniques on sperm parameters. *Syst Biol Reprod Med*. 2015;61(1):1-12.
387. Oktay K, Harvey BE, Partridge AH, Quinn GP, Reinecke J, Taylor HS et al. Fertility Preservation in Patients With Cancer: ASCO Clinical Practice Guideline Update. *J Clin Oncol*. 2018;36(19):1994-2001.
388. Corona G, Minhas S, Giwercman A, Bettocchi C, Dinkelman-Smit M, Dohle G et al. Sperm recovery and ICSI outcomes in men with non-obstructive azoospermia: a systematic review and meta-analysis. *Hum Reprod Update*. 2019;25(6):733-57.
389. Kuczyński W, Dhont M, Grygoruk C, Grochowski D, Wołczyński S, Szamatowicz M. The outcome of intracytoplasmic injection of fresh and cryopreserved ejaculated spermatozoa--a prospective randomized study. *Hum Reprod*. 2001;16(10):2109-13.
390. Nicopoullos JD, Gilling-Smith C, Ramsay JW. Does the cause of obstructive azoospermia affect the outcome of intracytoplasmic sperm injection: a meta-analysis. *BJU Int*. 2004;93(9):1282-6.

391. Hezavehei M, Sharafi M, Kouchesfahani HM, Henkel R, Agarwal A, Esmaeili V et al. Sperm cryopreservation: A review on current molecular cryobiology and advanced approaches. *Reprod Biomed Online*. 2018;37(3):327-39.
392. Li K, Rodriguez D, Gabrielsen JS, Centola GM, Tanrikut C. Sperm cryopreservation of transgender individuals: trends and findings in the past decade. *Andrology*. 2018;6(6):860-4.
393. Starmer BZ, Baird A, Lucky MA. Considerations in fertility preservation in cases of testicular trauma. *BJU Int*. 2018;121(3):466-71.
394. Kamischke A, Jürgens H, Hertle L, Berdel WE, Nieschlag E. Cryopreservation of Sperm From Adolescents and Adults With Malignancies. *J Androl*. 2004;25(4):586-92.
395. Bourne H, Richings N, Liu DY, Clarke GN, Harari O, Baker HW. Sperm preparation for intracytoplasmic injection: methods and relationship to fertilization results. *Reproduction, fertility, and development*. 1995;7(2):177-83.
396. Hawksworth DJ, Szafran AA, Jordan PW, Dobs AS, Herati AS. Infertility in Patients With Klinefelter Syndrome: Optimal Timing for Sperm and Testicular Tissue Cryopreservation. *Rev Urol*. 2018;20(2):56-62.
397. Kirkman-Brown JC, Martins MV. 'Genes versus children': if the goal is parenthood, are we using the optimal approach? *Hum Reprod*. 2020;35(1):5-11.
398. Parmegiani L, Accorsi A, Bernardi S, Arnone A, Cognigni GE, Filicori M. A reliable procedure for decontamination before thawing of human specimens cryostored in liquid nitrogen: three washes with sterile liquid nitrogen (SLN2). *Fertil Steril*. 2012;98(4):870-5.
399. Tedder RS, Zuckerman MA, Goldstone AH, Hawkins AE, Fielding A, Briggs EM et al. Hepatitis B transmission from contaminated cryopreservation tank. *Lancet*. 1995;346(8968):137-40.
400. Tomlinson M. Managing risk associated with cryopreservation. *Hum Reprod*. 2005;20(7):1751-6.
401. Perez-Sanchez F, Cooper TG, Yeung CH, al. e. Improvement in quality of cryopreserved human spermatozoa by swim-up before freezing. *Int J Androl*. 1994;17(3):115-20.
402. Machen GL, Harris SE, Bird ET, Brown ML, Ingalsbe DA, East MM et al. Utilization of cryopreserved sperm cells based on the indication for storage. *Investig Clin Urol*. 2018;59(3):177-81.
403. Liu S, Li F. Cryopreservation of single-sperm: where are we today? *Reprod Biol Endocrinol*. 2020;18(1):41.

404. Isachenko V, Rahimi G, Mallmann P, Sanchez R, Isachenko E. Technologies of cryoprotectant-free vitrification of human spermatozoa: asepticity as criterion of effectiveness. *Andrology*. 2017;5(6):1055-63.
405. Berkovitz A, Miller N, Silberman M, Belenky M, Itsykson P. A novel solution for freezing small numbers of spermatozoa using a sperm vitrification device. *Hum Reprod*. 2018;33(11):1975-83.
406. Li YX, Zhou L, Lv MQ, Ge P, Liu YC, Zhou DX. Vitrification and conventional freezing methods in sperm cryopreservation: A systematic review and meta-analysis. *Eur J Obstet Gynecol Reprod Biol*. 2019;233:84-92.
407. Sánchez, R., Risopatrón, J., Schulz, M., Villegas, J. V., Isachenko, V., & Isachenko, E. Vitrified sperm banks: the new aseptic technique for human spermatozoa allows cryopreservation at -86 °C. *Andrologia*. 2012; 44(6):433-435.
408. Schulz, M., Risopatrón, J., Uribe, P., Isachenko, E., Isachenko, V., & Sánchez, R. Human sperm vitrification: A scientific report. *Andrology*. 2020;8(6):1642-1650.
409. Cooper TG, Atkinson AD, Nieschlag E. Experience with external quality control in spermatology. *Hum Reprod*. 1999;14(3):765-9.
410. Cooper TG, Björndahl L, Vreeburg J, Nieschlag E. Semen analysis and external quality control schemes for semen analysis need global standardization. *Int J Androl*. 2002;25(5):306-11.
411. Matson PL. External quality assessment for semen analysis and sperm antibody detection: results of a pilot scheme. *Hum Reprod*. 1995;10(3):620-5.
412. Mortimer D, Björndahl L, Barratt CLR, Castilla JA, Menkveld R, Kvist U et al. A practical guide to basic laboratory andrology. 1st ed. Cambridge, UK ; New York: Cambridge University Press; In Press.
413. Neuwinger J, Behre HM, Nieschlag E. External quality control in the andrology laboratory: an experimental multicenter trial. *Fertil Steril*. 1990;54(2):308-14.
414. Mortimer ST, Mortimer D. Quality and Risk Management in the IVF Laboratory. 2nd ed. Cambridge: Cambridge University Press; 2015:402.
415. Westgard JO. Foreword to the second edition. In: Westgard JO, editor. Basic QC practices: training in statistical quality control for healthcare laboratories. Madison, WI, USA: QC Publishing; 2002.
416. Bohring C, Krause W. The intra- and inter-assay variation of the indirect mixed antiglobulin reaction test: is a quality control suitable? *Hum Reprod*. 1999;14(7):1802-5.
417. Kuster CE, Singer RS, Althouse GC. Determining sample size for the morphological assessment of sperm. *Theriogenology*. 2004;61(4):691-703.

418. Bland JM, Altman DG. Statistical methods for assessing agreement between two methods of clinical measurement. *Lancet*. 1986;1(8476):307-10.
419. Armitage P, Berry G, Matthews J. Statistical methods in medical research. Oxford: Blackwell Science; 2002:826
420. Motulsky H. Intuitive Biostatistics - A nonmathematical guide to statistical thinking. 4th ed. Oxford: Oxford University Press; 2017:608
421. Burtis CA, Ashwood ER. Tietz textbook of clinical chemistry. Philadelphia London Toronto Montreal Sydney Tokyo W.B. SAUNDERS COMPANY 1999:1917
422. Plaut DA, Westgard JOW. DC external quality assessment. In: Westgard JO, editor. Basic QC practices: training in statistical quality control for healthcare laboratories. Madison, WI, USA: QC Publishing; 2002. p. 125-63.
423. Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989.
424. Cooper TG, Noonan E, von Eckardstein S, Auger J, Baker HW, Behre HM et al. World Health Organization reference values for human semen characteristics. *Hum Reprod Update*. 2010;16(3):231-45.
425. Ford WC. Comments on the release of the 5th edition of the WHO Laboratory Manual for the Examination and Processing of Human Semen. *Asian J Androl*. 2010;12(1):59-63.
426. MacLeod J, Gold RZ. The male factor in fertility and infertility. II. Spermatozoon counts in 1000 men of known fertility and in 1000 cases of infertile marriage. *The Journal of urology*. 1951;66(3):436-49.
427. Boyd JC. Defining laboratory reference values and decision limits: populations, intervals, and interpretations. *Asian J Androl*. 2010;12(1):83-90.
428. Björndahl L, Barratt CL, Mortimer D, Jouannet P. 'How to count sperm properly': checklist for acceptability of studies based on human semen analysis. *Human reproduction (Oxford, England)*. 2016;31(2):227-32.
429. Ozarda Y, Sikaris K, Streichert T, Macri J, intervals ICoR, Decision L. Distinguishing reference intervals and clinical decision limits - A review by the IFCC Committee on Reference Intervals and Decision Limits. *Crit Rev Clin Lab Sci*. 2018;55(6):420-31.
430. Ceriotti F, Henny J. "Are my Laboratory Results Normal?" Considerations to be Made Concerning Reference Intervals and Decision Limits. *EJIFCC*. 2008;19(2):106-14.
431. Gomez E, Aitken J. Impact of in vitro fertilization culture media on peroxidative damage to human spermatozoa. *Fertil Steril*. 1996;65(4):880-2.

432. Mortimer D, Mortimer ST. Essential features in media development for spermatozoa, oocytes and embryos. In: Quinn P, editor. Culture Media, Solutions and Systems in Human ART. Cambridge: Cambridge University Press; 2014. p. 47-67.
433. Biggers JD, Whitten WK, Whittingham DG. The culture of mouse embryos in vitro. In: Daniel JC, editor. Methods in Mammalian Embryology. San Francisco, CA, USA: Freeman; 1971. p. 86-116.
434. Quinn P, Kerin JF, Warnes GM. Improved pregnancy rate in human in vitro fertilization with the use of a medium based on the composition of human tubal fluid. *Fertil Steril*. 1985;44(4):493-8.
435. Hotchkiss RS. Fertility in man. London: William Heineman Medical Books; 1945.
436. Eliasson R. Standards for investigation of human semen. *Andrologia*. 1971;3:49-64.
437. Lentner C. Geigy scientific tables. Vol. 1: Units of measurement, body fluids, composition of the body, nutrition. Basle: Ciba-Geigy Ltd; 1981:50.

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© 世界卫生组织 2021